


HDV 632

HERZOG
by **PAC**



The HDV 632 is a fully automated vacuum distillation Analyzer, which offers highest flexibility for the user. It can measure samples up to 650°C AET. Choose ASTM D1160 or define and save your own custom programs. Naturally, you can add, delete, and revise custom tests as your products and test needs evolve. The HDV 632 controls temperature, distillation

rate, and pressure with microprocessor precision. Volume measurements are made with digital accuracy, and test program selection is a simple “point and click” operation. Choose from an extensive library of built-in programs or your own predefined tests. The HDV 632 incorporates everything needed for fully automated distillation

testing, yet uses only 0.33 square meters (3.25 sq. ft.) of bench space. In fact, it even features built-in refrigeration systems — one for the cold trap plus a heater/cooler for the condenser bath — there is no need to run hoses to an external water supply or chiller system.

ADVANTAGES

- Fully automated D1160 vacuum distillation
- Unparalleled precision and accuracy
- Simple to operate with powerful HLIS Software
- Automatic cleaning procedure
- Automatic sample degassing and programmable vacuum ramp
- High safety standards

APPLICATION RANGE

- Volatility
- Liquid petroleum products
- Vacuum gasoil
- Biodiesel

STANDARD METHODS

- In accordance with:
 - ASTM D1160



solidpartners proven solutions

U.S.A.

PAC, LP | 8824 Fallbrook Drive | Houston, Texas 77064
T: +1 800.444.TEST | O: +1 281.940.1803 | F: +1 281.580.0719
sales.usa@pacpl.com | service.usa@pacpl.com

FRANCE

BP 70285 | Verson | 14653 CARPIQUET Cedex
T: +33 231 264 300 | F: +33 321 266 293
sales.france@pacpl.com | service.france@pacpl.com

GERMANY

Badstrasse 3-5 | P.O.Box 1241 | D-97912 Lauda-Königshofen,
T: +49 9343 6400 | F: +49 9343 640 101
sales.germany@pacpl.com | service.germany@pacpl.com

SINGAPORE

61 Science Park Road | #03-09/10 The Galen
Singapore Science Park III | Singapore 117525
T: +65 6412 0890 | F: +65 6412 0899
sales.singapore@pacpl.com | service.singapore@pacpl.com

NETHERLANDS

P.O.Box 10.054 | 3004 AB Rotterdam
Innsbruckweg 35 | 3047 AG Rotterdam
T: +31 10 462 4811 | F: +31 10 462 6330
sales.netherlands@pacpl.com | service.netherlands@pacpl.com

RUSSIA

Shabolovka Street | 34, Bldg. 2 | 115419 Moscow
T: +7 495 617 10 86 | F: +7 495 913 97 65
sales.russia@pacpl.com | service.russia@pacpl.com

CHINA

Room 1003, Sunjoy Mansion | No. 6 RiTan Rd.
Chao Yang District | Beijing 100020
T: +86 10 650 72236 | F: +86 10 650 72454
sales.china@pacpl.com | service.china@pacpl.com

INDIA

1508 | Dev Corpora | Pokhran Road No.1
Eastern Express Highway | Thane (W) - 400 601
T: +91-22-6700 4848 | F: +91-22-4228 4950
sales.india@pacpl.com | service.india@pacpl.com

MIDDLE EAST

A1 Quds Street, A1 Tawar road | LIU#H13 Dubai Airport Freezone
Near Dubai Airport (terminal 2) | P.O.Box #54781 | Dubai, UAE
T: +971 04 2947 995 | F: +971 04 2395 465
sales.middleeast@pacpl.com | service.middleeast@pacpl.com

SOUTH KOREA

#621 World Vision Building | 24-2, Youido-dong
Seoul 150-010
T: +82 2785 3900 | F: +82 2785 3977
sales.southkorea@pacpl.com | service.southkorea@pacpl.com

THAILAND

26th Floor, M. Thai Tower | All Seasons Place
87 Wireless Road | Lumpini, Phatumwan | Bangkok 10330
T: +66 2627 9410 | F: +662627 9401
sales.thailand@pacpl.com | service.thailand@pacpl.com

PAC Authorized Representatives are also located in most countries worldwide. For more information visit www.pacpl.com

HERZOG BY PAC

Herzog, originally established in 1937, is a long-established comprehensive line of laboratory instruments for testing distillation, flash point, vapor pressure, bitumen testing, cold flow properties, viscosity and other physical properties of materials.

HDV 632: HIGHLY PRECISE, RELIABLE ANALYSIS

EASY TO USE, PRECISION ANALYSIS

- Initiate tests with only a few clicks and keystrokes
- Digital Volume Measuring System with automatic calibration
- Dual transducer sensing system provides accurate, precise pressure measurements over a broad range of distillation pressures, including atmospheric pressures
- Digital, proportional pressure control eliminates transducer calibration and provides precise pressure regulation
- Easily removable receiver graduate eliminates need for any fluid connection
- Automatic cleaning sequence runs sufficiently long to clean the column and cool system for next distillation

ULTIMATE FLEXIBILITY, PERFORMANCE

- Custom distillation curves record results at additional volume and temperature points or stop testing at a specific volume/temperature point
- Allows pressure transducer placement at top of condenser column or at top of vapor column
- Enables receiving compartment temperature control at setpoints from +10° to +80°C
- Automatic sensor prevents sample from foaming into column during pressure reduction sequence
- Automatically calculates required sample weight based on sample density, and adjusts distillation results for weighing errors
- Provides results in tabular, graphic, or customized format
- Stores data internally or exports to LIMS

SPECIFICATIONS

Ordering Information	
	Model HDV 632 Automatic Vacuum Distillation Analyzer
Standard Test Method	
	ASTM D1160
Performance	
Sample Temperature Measurement	Auto calibrating sample thermometer circuit with integrated precision fixed resistors; 10 point on-screen calibration
Condenser Temperature Control:	0° to 80°C operating temperature; liquid jacket, independent of receiver; programmable for temperature ramping during distillations
Mechanized Volume Measurement	Digital infrared meniscus follower, unaffected by ambient light; precision stepper motor, resolution to 0.05 ml; calculated as function of motor steps, confirmed by optical detector; self-calibrating
Distillation Rate Control	Automatic—proven algorithm controls temperature in the flask; heater conforms to CE safety specifications; "Follow the heat curve" software function permits manual adjustment; distillate receiver provides independent temperature control with programmable temperature ramping during distillations
Documentation	
	Local display; parallel output port standard; RS-232 serial output port standard; PC and printer optional; compatible with HLIS for Windows®
Measurements	
Temperature	Non-inertial, low mass thermocouples protected by rigid metal thermowell for reliable operation
Pressure	Built-in gauges: differential (in-flask pressure) & barometric
Heating System	Percent recovered or evaporated calculated against recorded pressure variation in the flask during distillation run
Diagnostics & calibration	
	Real time status display and control for all mechanical and electrical systems locally or with optional PC network; 10 point on-screen calibration of resistance thermometer; recovery volume calibration with dry reference graduate
Physical	
Electrical Requirements	230 VAC (115 VAC with transformer), 3100 watts, 50/60 Hz
Dimensions	49 cm W x 63 cm D x 103 cm H (19.3 x 24.8 x 40.5 inches)
Weight	104 kg (230 lbs)

Copyright 2012/1 PAC L.P. All rights reserved 820-149

www.pacpl.com